

VETASSURE 18 L

Sterilisation Table

VETASSURE B18		VETASSURE S18				
DISPLAYED #, TEMPERATURE AND NAME	TYPE	TOTAL TIME	TYPE	TOTAL TIME	STERIZABLE LOAD Load in (.) is sterilizable only in B-Class	WRAPPED/ OPEN
1 HELIX - B&D TEST	TEST	25'	TEST	28' B&D only	N/A	✗
2 VACUUM TEST		15'		15'		✗
3 121 °C HOLLOW WRAPPED	B	40'	S	43'	Solid Small porous Simple hollow + (Narrow lumen) (Multiple casing)	✓
4 134 °C HOLLOW WRAPPED		34'		37'		✓
5 121 °C POROUS		47'		53'	Small sized porous + (Full porous)	✓
6 134 °C POROUS		40'		52'		✓
7 134 °C PRION	S	44'		47'	Load related to Creutzfeldt-Jakob disease	✓
8 134 °C RAPID		42'		33'	Solid Simple hollow + (Small porous)	✗

WEIGHT	Solids wrapped: 2.5 kg unwrapped: 4 kg	Solids wrapped: 1.5 kg unwrapped: 2.5 kg
	Porous: 1 kg	Porous: 0.5 kg

NOTE: Narrow lumen is an item open on one side with length equal to 5 times (or more) the diameter. Times are calculated in conditions of pre-warmed chamber.

TECHNICAL DATA

VETASSURE B18	VETASSURE S18
fractionated	pre and post vacuum
✓	
optional	
505Wx615Dx400H mm	
ø236x381.5 mm	
nr. 4 aluminium trays: 185Wx285Dx18H mm	
230 V ± 10%	
50-60 Hz	
2400 W	
OVERALL MEASURES	
INTEGRATED WATER PURIFICATION SYSTEM "ROSI"	
INTERNAL USB DATA LOGGER	
CONNECTION TO EXTERNAL DATA LOGGER	
WATER QUALITY SENSOR	
MTS MEMORY TEST SYSTEM	
BUILT-IN PRINTER	
INTEGRATED STEAM GENERATOR	
TYPE OF VACUUM	



VETASSURE B18