

VETASSURE 23 L

Sterilisation Table

VETASSURE B23		VETASSURE S23				
DISPLAYED #, TEMPERATURE AND NAME	TYPE	TOTAL TIME	TYPE	TOTAL TIME	STERIZABLE LOAD Load in (.) is sterilizable only in B-Class	WRAPPED/ OPEN
1 HELIX - B&D TEST	TEST	35'	TEST	28' B&D only	N/A	✗
2 VACUUM TEST		15'		15'		✗
3 121 °C HOLLOW WRAPPED	B	50'	S	53'	Solid Small porous Simple hollow + (Narrow lumen) (Multiple casing)	✓
4 134 °C HOLLOW WRAPPED		44'		47'		✓
5 121 °C POROUS		57'		63'	Small sized porous + (Full porous)	✓
6 134 °C POROUS		50'		62'		✓
7 134 °C PRION	S	54'		57'	Load related to Creutzfeldt-Jakob disease	✓
8 134 °C RAPID		42'		43'	Solid Simple hollow + (Small porous)	✗

WEIGHT	Solids wrapped: 3.5 kg unwrapped: 4.5 kg	Solids wrapped: 2.5 kg unwrapped: 2.5 kg
	Porous: 1.5 kg	Porous: 1 kg

NOTE: Narrow lumen is an item open on one side with length equal to 5 times (or more) the diameter. Times are calculated in conditions of pre-warmed chamber.

TECHNICAL DATA

VETASSURE B23	VETASSURE S23	
fractionated	pre and post vacuum	
		optional
		505Wx690Dx400H mm
		ø236x530 mm
		nr. 4 INOX trays; 185Wx440Dx19H mm
		230 V ± 10%
		50-60 Hz
		2400 W



VETASSURE B23